

Multichannel Electrochemical Measurement System

# HZ-Pro Series

- Electrode materials research
- Electrochemical analysis
- Surface processing research
- Corrosion protection research
- Metal materials research
- Catalyst research



HZ-Pro S4



HZ-Pro S12



HZ-Pro S2A

Introducing HZ-Pro, a modular, multichannel electrochemical measurement system engineered to support a wide range of measurements, from basic electrochemical measurements to corrosion tests, battery tests, and analytical tests.



## Features

- Capable of mounting a maximum of 12 potentiostat boards and booster boards (when using a 12-slot chassis)
- Offers an impedance measurement range up to 10MHz\*1
- Allows parallel connection of optional 5A booster boards (max. 55A)
- Performs simultaneous measurement of reference electrode, working electrode, and counter electrode potentials (simultaneous measurement of EIS)
- Supports external communication using digital input/output or PC communication functions (HoPath)
- Delivers automated temperature control with a thermostatic chamber (thermostatic chamber linkage function)

### Main unit

Lineup of 2/4/12-slot models

### Potentiostat board

Max. output of 12V, 500mA  
Supports external communication using digital input/output

### FRA


Provides EIS measurement (1MHz/10MHz\*1) with an additional license

### CPU board

Connects to a computer via 100BaseT Ethernet

### Booster board

Max. output of 10V, 5A  
Delivers up to 55A by connecting multiple boards




\*1 Requires a separately available potentiostat board with 10MHz FRA capability and a dedicated cell cable


# High-usability application software supports users with ample features

## Features

- Comprehensive coverage of measurement methods used in the electrochemistry, corrosion, and battery fields (p.10)
- Convenient electric quantity control (SOC, DOD, C-rate)
- Greater measurement flexibility (system configuration, additional conditions)
- Thorough safety features (power outage recovery, simplified measurements, etc.)
- “Channel Group function,” a powerful tool in the battery field
- User support features (language switching, automatic text saving, etc.)



Control of electrical current and voltage as the master channel  
Measurement of voltage by other channels  
Channel Group function




Control of multiple HZ-Pro units  
by multiple users  
HZ-Pro system configuration


# Data analysis supported by simple operations

## Features

- Various waveform processing functions (moving average, DC removal, etc.)
- Comprehensive range of analyses, including dV series and pause extraction series (GITT)
- Analysis methods defined by JIS, including internal resistance calculation
- Data search from sample information
- In-situ AC impedance analysis function (EIS software)
- Outputs in ZView format files



dQ/dV and dV/dQ curves (LiB 18650)



In-situ AC impedance analysis function

# HZ-Pro Product Lineup

## HZ-Pro chassis

HZ-Pro comes in three types of chassis with varying slot capacities to choose from to meet various applications and expansion needs. An optional external booster is also available.



2-slot chassis



4-slot chassis



12-slot chassis




Optional external booster

## Potentiostat boards and cell cables

The number of measurement channels can be expanded by adding potentiostat boards. The potentiostat board can be given FRA functionality (1MHz/10MHz<sup>\*1</sup>) with an additional license.



Potentiostat board (1MHz)



Cell cable (1MHz)



Potentiostat board (10MHz<sup>\*1</sup>)




Cell cable (10MHz<sup>\*1</sup>)

<sup>\*1</sup> Requires a separately available potentiostat board with 10MHz FRA capability and a dedicated cell cable.

## Booster board


The system's specifications can be expanded by adding booster boards. Booster boards can be connected in parallel, up to a maximum current of 55A.



Current booster board / Voltage booster board

## External system linkage function

The system can be linked to external devices via external analog and digital input/output.



## Technical Note

### FRA\*1 capable of measuring up to 10MHz for all-solid-state batteries

The potentiostat capable of electrochemical impedance measurement can measure impedance from 10  $\mu$ Hz to 10 MHz, the highest frequency range in the industry.





Image of an all-solid-state battery electrode and a Nyquist plot example

\*1 Requires a separately available potentiostat board with 10MHz FRA capability and a dedicated cell cable.



# Hoktnet Station

## Features

- Supports everything from creating measurement conditions to analysis and management
- Provides high usability
- Supports users with ample features
- Easy language switching operation (supports Japanese and English)

## Software configuration



### Hoktnet Client

Performs everything from creating measurement conditions to providing measurement control and generating measurement data.



### Hoktnet Analyzer

Displays measurement data in graph form and performs data analysis.




### Hoktnet Dat

Converts measurement data into CSV files.



### Hoktnet Explorer


Collectively manages measurement conditions, measurement data, and CSV files.





# Hoktnet Client

Allows flexible arrangements of measurement conditions with intuitive operation. Schematic diagrams visualize control waveforms, and a host of support features allow for detailed settings and customer-specific operating conditions.



# Hoktnet Dat

Converts measurement data into CSV files. The CSV file can be saved in a number of save formats, including Time Series, Cycle Series, Step Series, ZView Series, and EIS Series.



# Hoktnet Explorer


Collectively manages measurement conditions, measurement data, and CSV files. Desired files are found quickly by using advanced filtering and search functions, and past measurement data can be referenced from measurement histories. Important measurement data can also be organized together using the filing function.



# Hoktnet Analyzer

## Plot type


The data display can be switched with one click. Original plots can be created and registered among the measurement patterns.



## Bookmark




Point data is displayed with labels. By customizing the labels, information that cannot be displayed in 2D can be verified easily.

The labels can be registered in the point data.




## Multiple graphs

A multiple number of graphs (max. 9 graphs) can be displayed simultaneously. The linkage of point data allows easy data comparisons.




## Data series

Compatible with various data series, including dV series.



dQ/dV curves (LiB 18650)




dV/dQ curve (LiB 18650)

## Analysis function

Supports parameter designated analysis to intuitive analysis based on graph manipulation. Field-specific analyses can be easily executed.

## Sample making

Displays out-of-range data as sample information in a graph for easy identification of abnormal data.




No.	>電位(WE-RE 電位)	Y1軸(電流)
	mV	uA
× 54	529.540	49.714
× 55	539.459	49.698
× 56	549.368	49.693
× 57	559.602	49.680
× 58	569.524	49.698
× 59	579.438	49.662
× 60	589.358	49.705
× 61	599.586	49.696
× 62	609.509	49.724
× 63	619.436	47.169
× 64	629.350	27.867



# EIS

## Equivalent circuit and fitting

An equivalent circuit can be created and the value of each element determined with a simple operation.



# HZ-Pro Product Specifications

## Base unit


Item	Specification		
Model	HZ-ProS2A	HZ-ProS4	HZ-ProS12
No. of channel slots	2	4	12
External dimensions*1 W x H x D (mm)	370 x 150 x 480	256 x 350 x 500	440 x 350 x 500
Max. power consumption	250VA	600VA	1700VA
Weight*2	< 5.8kg	< 9.5kg	< 14.3kg
Power supply	Φ 1 AC100V~240V ±10% 50/60Hz		

\*1 Does not include protruding parts \*2 When no board is mounted


## PGS board

Item	Specification	
Model	HZA-PGS111	HZA-PGS101A
Max. output voltage	±12V	
Max. output current	±500mA	
Control voltage	±10V	
Min. control voltage resolution	160nV	
Current control resolution	0.004% of the range	
Voltage control accuracy	±0.03% of the setting ±1mV	
Current control accuracy	500mA~5μA range: ±0.03% of the setting ±0.08% of the range 500nA range: ±0.03% of the setting ±0.1% of the range 50nA range: ±0.03% of the setting ±0.16% of the range	
Voltage control response speed*3	2μS	
Voltage detection range	±10V, ±2.5V, AUTO	
Voltage detection accuracy	±0.03% of the reading ±1mV	
Input bias current	< 10pA	
Input impedance	> 1TΩ	
Current detection range	±500mA~±50nA (8 ranges), AUTO	
Current detection accuracy	500mA~5μA range: ±0.03% of the reading ±0.08% of the range 500nA range: ±0.03% of the reading ±0.1% of the range 50nA range: ±0.03% of the reading ±0.16% of the range	
Impedance (optional) Measurement frequency range	10μHz~10MHz	10μHz~1MHz
AUX	External analog control input, trigger IO × 2, analog input × 2, analog record output (V, I), stop input, hardware error output	
Weight	< 0.7kg	

\*3 At no load



HZA-PGS111 impedance accuracy map




HZA-PGS101A impedance accuracy map

## Booster

Item	Specification			
Model	HZA1005	HZA5001	HZAP0312	HZAP0612
Max. output voltage	±10V	±50V	±3V	-0.1V~+6V
Max output current	±5A	±1A	±100A	
Current range	±5A*1、±500mA	±1A*1、±100mA	±100A、±10A、±1A	
Current control accuracy	±0.2% of the range	±0.3%±500pA of the range	±0.3% of the range	
Current detection accuracy	±0.2% of the range	±0.3%±500pA of the range	±0.3% of the range	
Current control response speed	≒300μs*2	<50μs(No load)	≒600μs*3	
External dimensions*4 W x H x D (mm)	Inserted in the base unit	Inserted In the base unit	440 x 350 x 500	
Max. power consumption	130VA per board	115VA per board	1800VA	
Weight	<1.2kg	<1.2kg	<30kg	
Power supply	Supplied from the base unit	Supplied from the base unit	Φ1 AC 200V~240V±10% 50/60Hz	

\*1 Range x Number of parallels = Actual current range \*2 At a control value of 5A and load of 100mΩ \*3 At a control value of 100A and load of 1mΩ \*4 Excluding protruding parts



HZAP0312 and 0612 impedance accuracy map

## Measurement method

Item	Measurement method
Basic measurement	OCP (OCV), CA, CC, CP, CV, LSV, BE, IR, Electrochemical impedance spectroscopy Potential pulse, current pulse, potential function*4, current function*4
Charge/discharge measurement	CC mode, CCCV mode, CP mode, CCCP mode, GITT, CR mode (discharge only)
Corrosion measurement	REST potential, pitting potential, TAFEL, protective potential, zero shunt ammeter, reactivation rate, anodic polarization, crevice erosion repassivation potential
Analysis	NPV, DPV, SWV, LSSV*5, DPSV*5, SWSV*5
Other	Flat band potential, arbitrary waveform*6


\*4 Step and sweep measurements can be combined. \*5 Stripping voltammetry \*6 Potential and current control can be mixed. However, this accompanies a certain amount of step switching time.

## Analysis function

Item	Analysis method
Point data	Arbitrary points, coordinate specification, open circuit potential, measurement start potential, solution resistance
Battery & capacitor	DC internal resistance, short-circuit current (lead-acid battery), capacitor capacitance (linear approximation, power conversion), intersection internal resistance
CV analysis	Peak recognition, coulomb value, E1/2
Waveform analysis	Maximum/minimum value, average value, integration, length measurement, half-wave potential
Corrosion analysis	Pitting potential, reactivation rate
Physical properties evaluation	Flat band potential
Other functions	Auxiliary line drawing (between two points, linear approximation, horizontal line, vertical line, bisector), waveform processing (moving average, DC removal, difference, translation), data dump, bookmarks, user-defined plot types, zoom display, etc.

# Technical Note

## In-situ AC impedance analysis function




## Using HoPath\*1 to control HZ-Pro via a custom system \*1

Customers can control HZ-Pro from their own system. By setting the start and stop of measurements at any desired timing, a custom system incorporating the HZ-Pro can be realized.

\*1 Optional feature

### System configuration example (fuel cell evaluation system)



Inquiries

## HD MEIDEN HOKUTO

### Head office & Tokyo Office

4-22-13, Himonya, Meguro-ku, Tokyo, 152-003 Japan  
Email [honsha@meiden-hokuto.co.jp](mailto:honsha@meiden-hokuto.co.jp)(Tokyo)

### Osaka office

1-1-1 Nishinagasu-cho, Amagasaki-shi, Hyogo, 660-0805 Japan  
Email [osaka@meiden-hokuto.co.jp](mailto:osaka@meiden-hokuto.co.jp)(Osaka)

HOME PAGE <https://www.meidensha.co.jp/hkt>

■ Specifications are subject to change due to improvements in functions and performance.